

GERMS? No Thanks!







Sanitizes in 10 seconds!



296000 Birex® Quat™ Ready-to-use Spray 24 ounces | 12 per case



295999
Birex® Quat™ Ready-to-use Gallon
1 gallon refill | 4 per case



295507Birex® Quat™ Wipes 6 per case

Use Wipes, Liquids, or Both without Worry!

Compatible together. Compatible with over 30 surfaces.

- 2 minute TB kill time
- Over 30 EPA pathogen kill claims
- Sanitizes in 10 seconds. Eliminates 99.9% of germs on hard, non-porous surfaces
- Can be used on over 30 dental surfaces including dental chairs and plastics

DENTAL DISINFECTANT WIPE COMPARISON

	Birex® Quat™ Wipes	CaviWipes™	CaviWipes 1™	Super Sani- Cloth® Wipes	Optim® 33TB Wipes	Discide® Ultra Towelettes
Efficacy						
Benchmark Organisms Kill Time Including TB (min)	2	3	1	2	5	1
Sanitation (seconds)	10	-	-	-	30	-
Microorganisms Kill Claims	259	13	24	30	28	24
Product Specs						
Standard Wipe Size	6" x 6.75"	6" x 6.75"	6" x 6.75"	6" x 6.75"	6" x 7"	6" x 6.75"
Standard Wipe Count	160	160	160	160	160	160
Liquids Available	Yes	Yes	Yes	Yes	Yes	Yes
Compatibility						
Alcohol Percentage	21%	17%	23%	55%	Not Reported	63%
EPA Label Approximated Use Sites Approved	259	124	251	171	114	78
CoPolymers, Plastics, & Rubbers noted to be affected by exposure	-	Neoprene Discoloration	Neoprene Discoloration, Kraton & Incompatibility	Clear Plastic	-	Use care on Plexi-Glass
Metals noted to be affected by exposure	-	Brass, Aluminum, Carbon Steel & Copper	Brass & Copper Discoloration	Brass	Avoid prolonged exposure to Copper, Brass, Anodized Aluminum, some coatings	_
General Surfaces noted to be affected by exposure	-	_	Glass Spotting	Not intended for use on Touch Screen monitors, natural marble, unpainted wood	Not recommended on marble, polished wood, granite and leather	Use care on varnished surfaces

