Upgrade Your Surface Disinfectant

BirexSE is now BirexSE III: New and improved formula, same great benefits.



The Smart Way to Clean

- Spend up to 50% less with BirexSE III concentrate!
- Each packet makes 1 gallon of solution
- Superior material compatibility allows for high-frequency use on most surfaces*
- One-step cleaner effective even with 5% blood serum



BIPEXSE III EPA REG. NO 1043-128-51003					
PRODUCT	ITEM #	COUNT	SIZE		
NEW! BirexSE III Intro Pack	296041	6 packets	1 oz. p		

NEW! BirexSE III Intro Pack	296041	6 packets	1 oz. packets	128 oz.	6 gallons	
NEW! BirexSE III Operatory Pack	296042	12 packets	1 oz. packets	128 oz.	12 gallons	
NEW! BirexSE III Clinic Pack	296043	36 packets	1 oz. packets	128 oz.	36 gallons	
NEW! BirexSE III Accessory Pack	296040	2-32 oz. Bottles, 1-Gallon Mixing Jug, 1-Funnel, 2-Squirt Caps				

Note: BirexSE accessories should be discarded and replaced with BirexSE III accesory pack, as dilution ratios have changed.

NEW! Reusable Accessory Pack for seamless mixing and application

- 1 packet dilutes to 1 gallon
- Gallon jug and funnel for convenient preparation
- Diluted solution can be used for up to 35 days
- 2-32 oz. spray bottles for easy application
- Squirt cap to minimize aerosols and for disinfecting large areas

SAVE SPACE AND MONEY!

1 - 6ct BOX OF BirexSE III MAKES

6 GALLONS OF SURFACE DISINFECTANT



DILUTION

MAKES



MORE EFFICIENT DILUTION: 1 packet of BirexSE III = 1 gallon vs. 4 packets of BirexSE = 1 gallon

BirexSE III is effective against M. Bovis (TB), HIV-1, staphylococ-cus, salmonella, pseudomonas, streptococcus, and drug-resistant organisms such as MRSA and H1N1 (swine flu). BirexSE III is effective against Covid-19 (SARS-CoV-2) when used in accordance with the direction for use against Adenovirus Type 2. EPA List N, List B approved. BirexSE III complies with OSHA's Blood Borne Pathogens Standard and the CDC guidelines for dental disinfection and cleaning.

† Pricing Data from HenrySchein.com accessed April 28, 2021 **When compared to leading Ready-to-Use Surface Disinfectant Sprays.
*Can be used on most metal surfaces including chrome and stainless steel; appropriate for Vinyl, washable paints, and other sealed surfaces

